Sample Size Determination for a Flow Cytometry Assay

Elaine McVey, Friedrich Hahn, Megan Gottlieb, Ruiling Yuan, Perry Haaland

BD Technologies, Research Triangle Park, NC

ABSTRACT

Flow cytometry is an important assay technology that has been limited by a lack of computational tools appropriate for use by statisticians and bioinformaticians. In this application, we demonstrate the use of a new open source software package to perform standard statistical analyses on flow cytometry data. In this case, cell numbers available for assay were limited, and sample size planning was required. Data was available from a calibration experiment which mimicked the expected experimental conditions and provided large numbers of observations. Using the R package flowCore, part of the Bioconductor project. these flow cytometry data were resampled to obtain proper statistical estimates of assay sensitivity parameters for small sample sizes.

EXPERIMENTAL DATA

- A protein was identified as a surrogate marker for a phenotype of interest to be used in high throughput screening. Due to limits on time and cell availability, it was desirable to minimize sample size.
- Analysis was based on data from an artificial system using a mixture of two cells types. One cell was positive for the protein of interest and the other was negative. Mixtures were prepared representing 0.1, 5, 25, and 100 percent positive. This data was used to create a standard curve for the assav relating true percent positive to observed percent positive.

METHODS

Flow Analysis:

- The R Bioconductor package *flowCore* was used to read in raw FCS files, sample data points, and apply gating. Flow data was displayed using the flowViz package (Figures 1 and 2).
- Gates were set using the 1% positive samples. Separate gates were determined for each of four cell donors. These gates were applied to all samples (0%-100%) to determine the observed percent positive.

METHODS

Curve Fitting:

• The R package *calib* was used to fit four parameter logistic models with variance modeled as a power of the mean.

$$\mathsf{FPL}_{\mathsf{POM}}: \mathsf{y} = \frac{\beta_1 - \beta_2}{1 + (\frac{\mathsf{x}}{\beta_1})^{\beta_4}} + \beta_2 + \mu^{\theta} \, \sigma \varepsilon$$

Curves were fit to the standard curves based on the sampled data and detection limits were estimated (Figure 3).

- Detection Limits:
 - The reliable detection limit (RDL) was used to measure the sensitivity of the assay.
 - The RDL represents the minimum percent of positive cells that could be reliably detected.
 - . The calib package calculates the RDL by searching for the x-axis value for which the lower v-axis confidence limit is equal to the upper y-axis confidence limit at zero.

Resampling:

 At each of five hypothetical sample sizes (100, 250, 500, 1000, 2000 events), data was resampling fifty times, resulting in fifty standard curves and fifty estimates of the detection limit. Non-parametric confidence intervals were calculated based on these data, and the resulting dependence of sensitivity on sample size was determined (Figure 4).

RESULTS





flowCore Code - Reading from FCS files*

library(flowCore) fSet <- read.flowPlate(path = datafolder, wellAnnotation = anno.adf)

flowCore/flowViz Code - Gating*

library(flowViz) library(lattice) ## filter out debris debFilt <- rectangleGate('FSC.A' = c(0.2, 0.8), 'SSC.A' = c(0.0, 8)) xvplot(SSC.A ~ FSC.A, fSet.nondebris, smooth = FALSE, filter = debFilt) fSet.nondebris <- Subset(fSet, debFilt)</pre> ## gate positive cells boundPos.mat <- cbind(c(0, 0, .7, .7, 0), c(1.1, .35, .65, 1.1, 1.1))

colnames(boundPos.mat) <- c('SSC-A', 'PE-A') gatePos <- polygonGate('SSC-A', 'PE-A', filterId = 'PosThres', boundaries</pre> boundPos.mat)

xyplot(PE.A ~ SSC.A, fSet.nondebris, smooth = FALSE, filter = gatePos)

flowCore Code - Sampling*

sample 100 random events set.sub <- fsApply(fSet, function(x) {</pre> samp <- sample(1:10000, size = 100) new('flowFrame',exprs = x[samp,]) use.exprs = TRUE, simplify = FALSE) fSet.sub <- as(set.sub, 'flowSet')

calib Code - Curve Fitting**

0.6

0.8

library(calib)

fit.perc <- calib.fit(Result.df\$percentPos. Result.df\$obsPos. m = 3. type = 'fpl.pom')</pre> plot(fit.perc) fit.perc@rdl

* some flowCore code may not function as shown until the next public flowCore release (1.2) ** calib is expected to be posted to CRAN in early Fall 2007



Figure 3. Fitted standard curve with confidence limits for values sampled at one sample size (calib)

CONCLUSIONS AND DISCUSSION

- At least 1000 observations should be collected in order to have confidence that 5% positive cells would be reliably detectable.
- The availability of automated gating would make it possible to incorporate variability due to gating.
- Due to the binomial nature of the data, the standard deviation does not increase consistently with the mean, but we believe the detection limits are minimally effected.

REFERENCES

 Regression and calibration with nonconstant error variance. Chemometrics and Intelligent Laboratory Systems, Volume 9, Issue 3, December 1990, pp. 231-248, Marie Davidian and Perry D. Haaland.

ACKNOWLEDGMENTS

- For assistance with code.
 - Errol Strain flowCore and flowViz
 - Bvron Ellis flowCore
 - Daniel Samarov calib
- Pat McCutchen layout

2000

Figure 4. Relationship between detection limit

and sample size as determined by resampling

(50 iterations). Confidence intervals (95%) were

determined non-parametrically. The dotted line represents the assay sensitivity required for

future experiments.

EAM/pwm - BDT, 7/2007



100 80

Sample Size Determination for a Flow Cytometry Assay

- Elaine McVey, Friedrich Hahn, Megan Gottlieb, Ruiling Yuan, Perry Haaland

BD Technologies, Research Triangle Park, NC

ABSTRACT

Flow cytometry is an important assay technology that has been limited by a lack of computational tools appropriate for use by statisticians and bioinformaticians. In this application, we demonstrate the use of a new open source software package to perform standard statistical analyses on flow cytometry data. In this case, cell numbers available for assay were limited, and sample size planning was required. Data was available from a calibration experiment which mimicked the expected experimental conditions and provided large numbers of observations. Using the R package flowCore, part of the Bioconductor project. these flow cytometry data were resampled to obtain proper statistical estimates of assay sensitivity parameters for small sample sizes.

EXPERIMENTAL DATA

- A protein was identified as a surrogate marker for a phenotype of interest to be used in high throughput screening. Due to limits on time and cell availability, it was desirable to minimize sample size.
- Analysis was based on data from an artificial system using a mixture of two cells types. One cell was positive for the protein of interest and the other was negative. Mixtures were prepared representing 0.1, 5, 25, and 100 percent positive. This data was used to create a standard curve for the assav relating true percent positive to observed percent positive.

METHODS

Flow Analysis:

- The R Bioconductor package *flowCore* was used to read in raw FCS files, sample data points, and apply gating. Flow data was displayed using the flowViz package (Figures 1 and 2).
- Gates were set using the 1% positive samples. Separate gates were determined for each of four cell donors. These gates were applied to all samples (0%-100%) to determine the observed percent positive.

METHODS

Curve Fitting:

 The R package calib was used to fit four parameter logistic models with variance modeled as a power of the mean.

$$FPL_{POM} : y = \frac{\beta_1 - \beta_2}{1 + (\frac{x}{\beta_2})^{\beta_4}} + \beta_2 + \mu^{\Theta} \sigma \varepsilon$$

Curves were fit to the standard curves based on the sampled data and detection limits were estimated (Figure 3).

- Detection Limits:
 - The reliable detection limit (RDL) was used to measure the sensitivity of the assay.
 - The RDL represents the minimum percent of positive cells that could be reliably detected.
 - The calib package calculates the RDL by searching for the x-axis value for which the lower v-axis confidence limit is equal to the upper y-axis confidence limit at zero.

Resampling:

 At each of five hypothetical sample sizes (100, 250, 500, 1000, 2000 events), data was resampling fifty times, resulting in fifty standard curves and fifty estimates of the detection limit. Non-parametric confidence intervals were calculated based on these data, and the resulting dependence of sensitivity on sample size was determined (Figure 4).

R CODE EXAMPLES

flowCore Code - Reading from FCS files*

library(flowCore) fSet <- read.flowPlate(path = datafolder, wellAnnotation=anno.adf

flowCore/flowViz Code - Gating*

library(flowViz) library(lattice) ## filter out debris debFilt <- rectangleGate('FSC.A' = c(0.2, 0.8), 'SSC.A' = c(0.0, 8)) xvplot(SSC.A ~ FSC.A, fSet.nondebris, smooth = FALSE, filter = debFilt) fSet.nondebris <- Subset(fSet, debFilt)</pre> ## gate positive cells boundPos.mat <- cbind(c(0, 0, .7, .7, 0), c(1.1, .35, .65, 1.1, 1.1)) colnames(boundPos.mat) <- c('SSC-A', 'PE-A')

gatePos <- polygonGate('SSC-A', 'PE-A', filterId = 'PosThres', boundaries</pre> boundPos.mat)

xyplot(PE.A ~ SSC.A, fSet.nondebris, smooth = FALSE, filter = gatePos)

flowCore Code - Sampling*

sample 100 random events set.sub <- fsApply(fSet, function(x) {</pre> samp <- sample(1:10000, size = 100) new('flowFrame',exprs = x[samp,]) use.exprs = TRUE, simplify = FALSE) fSet.sub <- as(set.sub, 'flowSet')

calib Code - Curve Fitting**

library(calib)

fit.perc <- calib.fit(Result.df\$percentPos. Result.df\$obsPos. m = 3. type = 'fpl.pom')</pre> plot(fit.perc) fit.perc@rdl

* some flowCore code may not function as shown until the next public flowCore release (1.2) ** calib is expected to be posted to CRAN in early Fall 2007





dence limits for values sampled at one sample size (calib)

CONCLUSIONS AND DISCUSSION

- At least 1000 observations should be collected in order to have confidence that 5% positive cells would be reliably detectable.
- The availability of automated gating would make it possible to incorporate variability due to gating.
- Due to the binomial nature of the data, the standard deviation does not increase consistently with the mean, but we believe the detection limits are minimally effected.

REFERENCES

 Regression and calibration with nonconstant error variance. Chemometrics and Intelligent Laboratory Systems, Volume 9, Issue 3, December 1990, pp. 231-248, Marie Davidian and Perry D. Haaland.

ACKNOWLEDGMENTS

- Errol Strain flowCore/flowViz
- Byron Ellis flowCore
- Daniel Samarov calib
- Pat McCutchen lavout

JBA/owm - BDT, 1/2007



Figure 3. Fitted standard curve with confi-



Sample Size (pre-gating)

Figure 4. Relationship between detection limit and sample size as determined by resampling (50 iterations). Confidence intervals (95%) were determined non-parametrically. The dotted line represents the assay sensitivity required for future experiments.